medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43-512 57 95 15
F +43-512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium bvba
Posthoornstraat 13/1
3582 Koersel
Belgium
T: + 32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Brasil
Rua Itapeva 538 cj. 42
Bela Vista, São Paulo - SP
CEP: 01332-000
Brasil
T +55 (11) 3500 8000
sac@medibrasil.com
www.medibrasil.com

medi Danmark ApS
Vejlegardsvej 59
2665 Vallensbaek Strand
Denmark
T +45 70 25 56 10
F +45 70 25 56 20
kundeservice@sw.dk
www.medidanmark.dk

medi Bayreuth España SL
C/Canigó 2 – 6 bajos
Hospitalet de Llobregat
08901 Barcelona
Spain
T +34 93 260 04 00
F +34 93 260 23 14
medi@mediespana.com
www.mediespana.com

medi France
Z.I. Charles de Gaulle
25, rue Henri Farman
93297 Tremblay en France Cedex
France
T +33 1 48 61 76 10
F +33 1 49 63 33 05
infos@medi-france.com
www.medi-france.com

medi Hungary Kft.
Bokor u. 21.
1037 Budapest
Hungary
T +36 1371 0090
F +36 1371 0091
info@medi.hu
www.medi.hu

medi Nederland BV
Heusing 5
4817 ZB Breda
The Netherlands
T +31 76 57 22 555
F +31 76 57 22 565
info@medi.nl
www.medi.nl

medi Orient M.E.A. ZZC
TAURUS GROUP
Suite No. M01 Wasel Al Mamzar Bldg.
Al Wuhida Road, Al Mamzar
P.O.BOX NO. 91464
DUBAI, UAE
United Arab Emirates
T: +971-4 2557344
F: +971- 4 2557399
md@taurusgroup.net

medi Polska Sp. z.o.o.
Zygmunta Starego 26
44-100 Gliwice
Poland
T +48-32 230 60 21
F +48-32 202 87 56
info@medi-polska.pl
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
Rua do Centro Cultural, no. 43
1700-106 Lisbon
Portugal
T +351 21 843 71 60
F +351 21 847 08 33
medi.portugal@medibayreuth.pt
www.medi.pt

4 051526 890664

E013521/10.2017

medi Epibrace

Gebrauchsanweisung. Instructions for use. Mode d'emploi. Instrucciones de uso. Instruções para aplicação. Istruzioni per l'uso. Gebruiksaanwijzing. Brugsvejledning. Bruksanvisning. Návod k použití. Uputa za uporabu. Инструкция по использованию. Kullanma kılavuzu. Instrukcja zakładania. Οδηγία εφαρμογής. Használati útmutató. Інструкція з використання. Navodila za uporabo. Instrucţiuni de utilizare.

medi. I feel better.

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Important notes
This medical product is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability according to the Medical Devices Act will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the support, please consult your doctor or orthotist immediately. Do not wear the brace over open wounds, and use it only as instructed by your doctor or orthotist.

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medi Epibrace

Intended purpose
medi Epibrace is an epicondylitis brace. The product is exclusively to be used for the orthotic fitting of the muscle base and only on unbroken skin.

Indications
Excessive strain on the tendon tissue at the end of the elbow, i.e. Tennis elbow (Epicondylitis humeri radialis) or golfer's elbow (Epicondylitis humeri ulnaris)

Contraindications
No side effects have been reported to date. In the following circumstances you should consult your doctor before use:

- Disorders or injuries of the skin in the area where the product is used, above all in case of signs of inflammation (excessive heat, swelling or reddening).
- Sensory or circulatory disturbances of the forearm or the hand.
- Disorders of lymphatic drainage and swelling of unclear origin of the soft tissues distal to the elbow.
- Impingement syndrome (pressure irritation of the nerves supplying the arm).

If you have a known history of contact allergy or tenderness to pressure, you should wear a thin cotton armsleeve under the medi Epibrace as a preventive measure. In principle, supports or braces that fit too tightly to or around the body can cause local pressure phenomena or, more rarely, to constriction of blood vessels or nerves.

Wearing directions
During longer periods of rest (e.g. when sleeping) you should take the medi Epibrace off.

Applying the strap

- Open the strap on the brace and lie it down, fully unrolling it, in front of you. Position your forearm on the brace with your palm facing upwards so that the light grey positioning mark on the pad is positioned on the side that

you feel pain. Then attach the pad to the positioning mark so that the tip of the pad is parallel to the elbow joint (Fig. 1). Example: If you have tennis elbow (Epicondylitis humeri radialis) in your right arm, the positioning mark should now be on the right side of your right arm.

- Thread the strap through the plastic loop and attach the fastening. The distance between the edge of the fastening and the outer elbow should be approx. 3 cm when the elbow is bent (Fig. 2).
- Refer to the table for guidance values on how tight the fastening can be done up (Fig. 3). Your forearm circumference (measured 5 cm below the elbow joint) should be used as the reference value. Tighten the strap until you attain the value shown on the chart. You can adjust the pressure according to your wishes.
- The clasp should now be firmly positioned on your arm, without slipping or constricting the arm. Form a fist with your hand. You should be able to feel your muscle under the pad. When your muscles are relaxed, the brace should sit against the pad without a noticeable feeling of pressure.
- In principle, the product should be used in accordance with medical supervision and should be fitted by a trained specialist, i.e. an orthopaedic technician, the first time that you put it on. If you wish to use other products with this product, then this should be discussed and agreed with your GP.
- If your arm or hand becomes numb, loosen the aid or remove it if necessary. If your condition persists, consult a doctor.

Care instructions

Velcro fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions
Store in a dry place at room temperature. Protect from heat, direct sunlight and moisture.

Material composition
Polyester, polyamide, polyurethane, polypropylene

Liability
The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal
The product can be disposed of in the domestic waste.

Your medi team
wishes you a quick return to full fitness.

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